



MULTI-LIGAND CONTROL-TRI LEVEL PRODUCT CODE: ML-300
LOT#:MLAC1E1 EXP: 2024-05-03

INTENDED USE
The Multi-Ligand Controls are intended for use as an assayed quality control material to monitor the consistency of performance of laboratory test procedures associated with determination and monitoring of the clinical status.

SUMMARY AND EXPLANATION
The use of quality control material to assist in the assessment of precision in the clinical laboratory is an integral part of laboratory practices. Controls that contain varied levels of analytes are necessary to insure precision and accuracy in immunosay systems.

REAGENTS
Monobind's Multi-Ligand Controls are intended to be used in the exact manner as patient samples. The control is packaged as 6 vials of 3.0 ml, dried. The analyte activities are adjusted to concentrations in the low, middle and high range in order to monitor the efficacy of the procedure in use.

INSTRUCTIONS FOR USE
1) Bring the vials to room temperature before use.
2) Carefully unscrew and remove cap.
3) Add three (3) ml of distilled or deionized water to each vial. Close the cap tightly and let the contents mix thoroughly for 30 minutes
4) Aliquot the materials in 0.5 ml aliquots in cryo vials and store at -20°C.

STORAGE, STABILITY AND DISPOSAL
This product will be stable until the expiration date when stored unopened at 2 to 8°C. Once the control is reconstituted, all analyses will be stable for 7 days when stored tightly capped at 2 to 8°C with the following exceptions: 1) C-Peptide should be assayed immediately after reconstitution, and 2) Folate and Insulin will be stable for 1 day. To avoid contamination, it is recommended labs aliquot required quantities into vials before each use.

After reconstituting, controls should be tightly capped and returned to refrigerator 2 to 8° C as soon as practical after usage. (Long term room temperature storage is not supported.) After reconstituting, controls should be tightly capped and frozen within 2-hours. Once thawed, do not refreeze the control; discard remaining material. It is recommended that customers aliquot control into separate containers before freezing to allow for usage on different days. Outdated material should be discarded as a biohazardous component.

Table with 3 columns: STORAGE, STABILITY, TEMPERATURE. Rows include Lympholized Unopened, Reconstituted, Opened, and Reconstituted, Opened.

EXPECTED RANGE OF VALUES
The mean values printed in this insert were derived from replicate analyses and are specific for this lot of product. The tests listed were performed by Monobind QA using representative lots of this product, as well as those of Monobind's AccuBind® ELISA and AccuLite® CLIA reagents.

Individual laboratory means should fall within the corresponding acceptable range; however laboratory means may vary from the listed values during the life of this control. Therefore, each laboratory should establish its own means and acceptable ranges for the product used, using Monobind's assignment only as guide. A trend log should be maintained for batch to batch consistency of the test. Variations over time and between laboratories may be caused by a) differences in laboratory personnel, b) improper technique, c) instrumentation and reagents, d) improper dilutions from the stated manufacturer's procedure, and/ or e) modifications in the manufacturer's test procedure.

Refer to http://www.monobind.com/site/qc-documents.html for any updated insert information.

WARNING AND PRECAUTIONS
FOR IN VITRO DIAGNOSTIC USE

All products that contain human serum have been found to be non-reactive for HIV 1&2, HIV-Ag, HBsAg, HCV and RPR by FDA required tests. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-395.

Revision: 2 Date: 2023-03-08 Product Code: ML-300

For Orders and Inquiries, please contact

Monobind Inc.
100 North Pointe Drive
Lake Forest, CA 92630 USA

Tel: +1 949.951.2565 Mail: info@monobind.com
Fax: +1 949.951.3039



Capartners® EL, Eastman's 13
9961 (DM) Main, The Netherlands
www.capartners.eu

Please visit our website to learn more
about our products and services.

Glossary of Symbols
EN ISO 15189



In Vitro
Diagnostic
Medical
Devices



Control
Material
Manufacturer

Used By
(Equipment Only)

Authorized Resale
in European Country

DOCUMENT HISTORY
PREPARED BY: MWP DEPT: QC VERIFIED BY: Ashok DEPT: QA
APPROVED BY: [Signature] DEPT: Administration EFFECTIVE DATE: REVISION: 2
DCO:

Table with columns: Analyte, Range, Range, Range, Method. Rows include categories like Allergy, Anemia, Bone Metabolism, Cancer Markers, Cardiac Markers, Diabetes, Fertility, Growth Deficiency, Steroids, Thyroid, and Thyroid VAST.