



Semi-Quantitative ELISA Antibody Tests

Routine screening for TORCH Antibodies during pregnancy or to monitor vaccine- response with excellent assay performance:

TORCH ELISA	Sensitivity True Positives	Specificity True Negatives	Samples Tested
Toxo IgG	100%	99.1%	120
Toxo IgM	100%	96.1%	135
Rubella IgG	98.9%	100%	100
Rubella IgM	100%	97.4%	130
CMV IgG	95.2%	95.8%	131
CMV IgM	84.6%	84.1%	126
HSV 1&2 IgG	85.4%	86.7%	127
HSV 1&2 IgM		<i>Coming Soon</i>	

- 
PLASMA & SERUM
 Detect, monitor treatment and recovery, as well as possible immunity*
- 
HIGHLY ACCURATE
 Well-established cutoff (> 10 RV) using 120 patients traceable to international standards as available for
 - Toxoplasma (NIBSC Code: 13/132)
 - Rubella (NIBSC Code: 67/182 and RUBI-1-94)
- 
REPRODUCIBLE RESULTS
 Uses recombinant antigens or a single viral line (Rubella IgG) with monoclonal mu-capture or IgG specific antibodies
- 
EASY & CONVENIENT PROCEDURE
 Ready-to-use reagents, no shaking, and room-temperature incubation
- 
FAST RESULTS
 75-80 minutes

* According to Structure and Function of Immunoglobulins, [READ MORE](#)



Scan me!


sales@monobind.com

www.monobind.com




Additional Benefits



HIGH YIELD

One-point cut-off determination allows more patients samples for a lower cost-per-test



EASILY AUTOMATED

Adaptable to open systems, including those sold by Monobind**

** *Monobind offers a line of ELISA analyzers*
[VIEW INSTRUMENTS](#)

ITEM

- | | |
|------------------|--|
| 6125-300A | Anti-Toxoplasma (Toxo) IgG
AccuBind ELISA Kit - 96 wells |
| 6225-300A | Anti-Toxoplasma (Toxo) IgM
AccuBind ELISA Kit - 96 wells |
| 6325-300A | Anti-Rubella (Rubella) IgG
AccuBind ELISA Kit - 96 wells |
| 6425-300A | Anti-Rubella (Rubella) IgM
AccuBind ELISA Kit - 96 wells |
| 6525-300A | Anti-Cytomegalovirus (CMV) IgG
AccuBind ELISA Kit - 96 wells |
| 6625-300A | Anti-Cytomegalovirus (CMV) IgM
AccuBind ELISA Kit - 96 wells |
| 6725-300A | Anti-Herpes Simplex Virus 1&2
(HSV 1&2) IgG
AccuBind ELISA Kit - 96 wells |
| 6825-300A | Anti-Herpes Simplex Virus 1&2
(HSV 1&2) IgM
AccuBind ELISA Kit - 96 wells
(Coming Soon) |

Contact us today to learn more



www.monobind.com



sales@monobind.com



Scan me!

