

# Semi-Quantitative ELISA Antibody Tests

Routine screening for TORCH Antibodies during pregnancy or to monitor vaccine- response with excellent assay performance

#### PLASMA & SERUM

Detect, monitor treatment and recovery, as well as possible immunity\*

\* According to Structure and Function of Immunoalobulins

### HIGHLY ACCURATE

Well-established cutoff (→→ 10 RV) using 120 patients traceable to international standards as available for -Toxoplasma (NIBSC Code: 13/132) -Rubella (NIBSC Code: 67/182 and RUBI-1-94)

## REPRODUCIBLE RESULTS

Uses recombinant antigens or a single viral line (Rubella IgG) with monoclonal mu-capture or IgG specific antibodies

## EASY & CONVENIENT PROCEDURE

Ready-to-use reagents, no shaking, and roomtemperature incubations

### FAST RESULTS

75-80 minutes

Monobind Inc.

www.monobind.com





### **Additional Benefits**

### **HIGH YIELD**

√ Manual method requires only onepoint cut-off determination allows more patients samples at a lower cost-per-test.

#### **EASILY AUTOMATED**

√ Adaptable to open systems, including those sold by Monobind\*\*

\*\* Monobind offers a line of ELISA analyzers SEE WEBSITE FOR MORE





6125-300A	Anti-Toxoplasma (Toxo) IgG AccuBind ELISA Kit - 96 wells	6525-300A	Anti-Cytomegalovirus (CMV) IgG AccuBind ELISA Kit - 96 wells
6225-300A	Anti-Toxoplasma (Toxo) IgM AccuBind ELISA Kit - 96 wells	6625-300A	Anti-Cytomegalovirus (CMV) IgM AccuBind ELISA Kit - 96 wells
6325-300A	Anti-Rubella (Rubella) IgG AccuBind ELISA Kit - 96 wells	6725-300A	Anti-Herpes Simplex Virus 1&2 (HSV 1&2) IgG
6425-300A	Anti-Rubella (Rubella) IgM AccuBind ELISA Kit - 96 wells	6825-300A	AccuBind ELISA Kit - 96 wells Anti-Herpes Simplex Virus 1&2 (HS)

Anti-Herpes Simplex Virus 1&2 (HSV 1&2) IgM AccuBind ELISA Kit - 96 wells

(Coming Soon)



Photo for example only. Not actual kit.



Contact us today to learn more



